LIMITED WARRANTY

Your Medline product is warranted to be free of defects in materials and workmanship for the duration of use by the original purchaser. This device was built to precise standards and carefully inspected prior to shipment. This limited warranty is an expression of our confidence in the materials and workmanship of our products and assures years of convenient, dependable service. In the event of a defect covered by this warranty, we will, at our option, repair or replace the device. This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters and grips, which are subject to normal wear and replacement.

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Medline Canada
1-800-396-6996
www.medline.ca | Canada@medline.com

Medline Mexico
01-800-831-0898
www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information. We reserve the right to correct any errors that may occur within this manual. ©2010 Medline Industries, Inc. Guardian Signature is a trademark of Medline Industries, Inc.

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GUARDIAN
• S I G N A T U R E •

PADDED TRANSFER BENCH WITH COMMODO OPENING

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READ ALL INSTRUCTIONS BEFORE USING.

WARNING

When transferring to or from the bench, do not put your entire weight on the front edge of the seat. Failure to comply might result in the seat pad popping off.

WARNING

Please read carefully before use. Failure to comply may result in severe injury.

WARNING

Ensure all detent push buttons are properly engaged in position prior to use.

WARNING

Do not exceed the maximum weight capacity listed on your product.

WARNING

Ensure suction cups are securely affixed to the tub floor. The suction cup legs must be secured to a clean level surface.

WARNING

Before use check to make sure screws are tight and product is stable.

FAILURE TO FOLLOW INSTRUCTIONS MAY RESULT IN SERIOUS INJURIES.

SAVE THESE INSTRUCTIONS.

ASSEMBLY

Backrest

- Determine which direction the user will face when sitting on the bench in the tub.
- To attach the backrest in the appropriate direction, slide the backrest poles into the tube openings underneath the seat assembly.
- Depress the detent push buttons and engage completely into holes to lock in place.

Commode Seat (Model 98013A only)

- To change the direction of the commode seat, pull up on seat to release clips.
- Rotate the seat, align clips on frame and press down making sure it is securely attached. Ensure commode opening faces front.

INSTALLATION

- Place the transfer bench over the tub wall with the side arm and suction footpieces inside the tub.
- The inner legs of the transfer bench should be as close to the far wall of the tub as possible.
- Adjust legs to desired height and to make level by depressing detent buttons and moving to desired position.
- Test the bench to ensure it is level and all legs are securely attached.
- Push down gently on the transfer bench to securely affix the suction cup feet to the tub floor.
- Test the unit for stability by vigorously attempting to move the transfer bench.

USING YOUR TRANSFER BENCH WITHOUT ASSISTANCE

Be sure to turn slowly, without twisting, throughout the following procedures to prevent injuries.

1. Transferring from a seated position:

- Place your wheelchair next to bench, facing front of the bathtub and engage wheel locks.
- Transfer your torso to the transfer bench first, putting your body weight in the middle of the seat, not on the seat edge. Lift your legs into the tub. Be sure to turn slowly, without twisting, to prevent injuries.
- After bathing lift your legs out of the tub first, then your torso.

2. Transferring from a standing position:

- Sit in the middle of the seat panel on the side of the bench that is outside of the tub.
- Extend your arms to grab the side arm of the bench for support and slide your body over to the center of the chair. Lift each leg over the tub wall.
- Slide over to a comfortable position in center of tub.
- After bathing, move your legs out of the tub, turn slowly and firmly grasp the backrest and stand.

NOTE– Tips On This Product Contain Dry Natural Rubber.

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